

RealGUIDE™ Z3D Guided Surgery

Designed for ZimVie's TSX™, TSV®, Trabecular Metal™, and Eztetic® Dental Implants



FULLY
GUIDED



RealGUIDE™ Z3D GUIDED SURGERY KIT

The new RealGUIDE Z3D Guided Surgery Kit contains instruments needed for the planned surgery. It is designed to work with ZimVie's **TSX, TSV, Trabecular Metal, and Eztetic** Implant Systems.

The kit information is available in the RealGUIDE Software Suite, an open system for diagnosis and advanced surgical guide design that will support you with customized drilling reports to be well-prepared for dental implant therapy.



RealGUIDE Software Suite
with PLAN, APP, GUIDE Modules

- Fully guided placement of ZimVie Dental Eztetic 3.1 mm Implant. TSX 3.1 mm, 3.7 mm, 4.1 mm, 4.7 mm and TSV/TM Implants in sizes 3.7 mm, 4.1 mm, and 4.7 mm.
- Color-coded for ease-of-use.
- Compact kit, which is useful for space-constrained dental offices.
- No handles or keys are necessary while preparing the osteotomy.
- DLC (Diamond-like Carbon) coating on the drills provides wear resistance characteristics.
- Includes implant mounts hex feature to align with the hex feature on the corresponding guide sleeves.
- Staging pedestals for holding the implant vial to remove the Fixture Mount Transfer (FMT) of the TSV implants or to mount the TSX implants is available. This enables the TSV implants to be placed fully guided through the surgical guides.
- Both soft and dense bone protocols are available.

Product	Part No.
RealGUIDE Z3D Surgical Kit	3DM0070ZUS



MORE INFORMATION
ABOUT THE RealGUIDE Z3D
GUIDED SURGERY KIT AT
ZIMVIE.COM

For more information, visit ZimVie.com

ZimVie
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Phone: +1-561-776-6700
Fax: +1-561-776-1272


Restoring Daily Life.®

Unless otherwise indicated, as referenced herein, all trademarks are the property of ZimVie; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Holdings, Inc., and distributed and marketed by ZimVie Dental (and, in the case of distribution and marketing, its authorized marketing partners). The RealGUIDE Z3D is manufactured by 3Diemme, a subsidiary of ZimVie. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV1214 REV A 05/23 ©2023 ZimVie. All rights reserved.

