

Dríva™ Gold Series Drills



Designed for high performance and excellent cutting efficiency, Dríva Gold Series Drills are reusable, externally irrigated, and compatible with existing surgical protocols for TSX®, Trabecular Metal®, Tapered Screw-Vent®, and Eztetic® Implants.

1. Efficient

- Deep cutting flutes designed for improved cutting efficiency compared to the original Dríva Drills
- Biocompatible ceramic coating increases hardness and durability

3. Convenient

- Solid construction facilitates easy cleaning and sterilization
- Gold color provides excellent visibility of grooves and bands that indicate drilling depth

2. Compatible

- Familiar soft- and dense-bone protocol are used for osteotomy preparation
- Black axial stripe on select drills indicates compatibility with Drill Stops

Ordering Information

Driva Gold Series Drills for TSX, Trabecular Metal, and TSV Implants

Diameter	Drill Length	
	16 mmL	22 mmL
2.3 mmD	TSV23SG	TSV23G
2.8 mmD	TSV28SG	TSV28G
3.4 mmD	TSV34SG	TSV34G
3.8 mmD	TSV38SG	TSV38G
4.4 mmD	TSV44SG	TSV44G
5.1 mmD	TSV51SG	TSV51G

Driva Gold Series Step Drills for TSX, Trabecular Metal, and TSV Implants

Diameter	Drill Length	
	16 mmL	22 mmL
3.4/2.8 mmD	TSV34D28SG	TSV34D28G
3.8/3.4 mmD	TSV38D34SG	TSV38D34G
4.4/3.8 mmD	TSV44D38SG	TSV44D38G
5.1/4.4 mmD	TSV51D44SG	TSV51D44G
5.7/5.1 mmD	TSV57D51SG	TSV57D51G

Driva Gold Series Surgical Kits for TSX, Trabecular Metal, and TSV Implants

Description	Item No.
Tapered Screw-Vent Surgical Kit, Complete	TSVKITG
Tapered Screw-Vent Starter Kit (3.7 and 4.7 mmD only)	STRKITG
Surgical Tray and Driva Gold Series Drills for TSV only	TSVTRAYDG

*Available for a limited time only, while supplies last.


Driva Gold Series Module and Drills for Eztetic Implants

Description	Item No.
Eztetic Surgical Module, Driva Gold Series	NPMODG
Driva Gold Series Drill, 2.4 mmD, 22 mmL	EZT24G
Driva Gold Series Step Drill 2.8/2.4 mmD, 22 mmL	EZT28D24G



Enhanced Cutting Efficiency¹

Original vs. Gold Series
(2.3 mm Driva Drill)



¹Data on file, available on demand

For more information, visit ZimVie.com

ZimVie
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Phone: +1-561-776-6700
Fax: +1-561-776-1272
dentalCS@ZimVie.com

Biomet 3i (UK & Ireland) Ltd.
Reading Business Centre
Suite 807, 8th Floor Fountain House
2 Queens Walk, Reading, Berks
RG1 7QF, United Kingdom
Telephone (UK): 0800 652 1233
Telephone (Ireland): 1 800 552 752
ZV.UKorders@ZimVie.com

ZimVie
Restoring Daily Life.®