


Absorbable collagen
membranes

BioMend[®] and BioMend Extend[™]



1. Absorbable¹

Eliminates second-stage surgery for membrane removal, designed to reduce wound trauma and surgical chair time

2. Cell-Occulsive

Serves as barrier to prevent epithelial cell migration and allow passage of essential nutrients¹

3. Space Maintaining

A rigid enough scaffold for tissue regeneration in GTR and GBR procedures



BioMend Membrane being prepared for placement in the sinus cavity

 **ZimVie**

BioMend® and BioMend Extend™

The absorbable membrane rigid enough to create and maintain space

Versatile Barrier Function

- Indicated for use in guided tissue and guided bone regeneration procedure
- BioMend Membrane is designed to be absorbed in approximately 8 weeks¹
- BioMend Extend Membrane is designed to be absorbed in approximately 18 weeks¹

Convenient Handling

- Tear resistant, suturable, pliable
- Easy to handle even when hydrated
- Rigid enough to create and maintain space



1. Exposed buccal wall defect on #3.



2. Membrane being prepared for placement.



3. Membrane cut for placement.



4. Membrane placed in socket.



5. Sutures placed to close site.



6. Final aesthetic result.

Clinical photos ©2013 Gary T. Jones, DMD, Sandhills Oral and Maxillofacial, Fayetteville, NC, USA. All rights reserved. Individual results may vary.
¹ Data on file with Collagen Matrix, Inc.

Ordering Information

Catalog Number

Description

0103Z BioMend Membrane, 15 x 20 mm

0105Z BioMend Membrane, 20 x 30 mm

0107Z BioMend Membrane, 30 x 40 mm

Catalog Number

Description

0140Z BioMend Extend Membrane, 15 x 20 mm

0141Z BioMend Extend Membrane, 20 x 30 mm

0142Z BioMend Extend Membrane, 30 x 40 mm



For more information, visit ZimVie.com

ZimVie
 4555 Riverside Drive
 Palm Beach Gardens, FL 33410
 1-800-342-5454
 Phone: +1-561-776-6700
 Fax: +1-561-776-1272

ZimVie
 Restoring Daily Life.®

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. In certain countries, BioMend and BioMend Extend Membranes are manufactured by Collagen Matrix, Inc. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV1662 REV A 10/23 ©2023 ZimVie. All rights reserved.

