



Strong Meets
Supple

Copios[®] Pericardium Membrane



THE UNIQUE
TUTOPLAST[®]
PROCESS

Demonstrated Performance

- Clinically demonstrated performance in guided bone regeneration procedures,^{1,2} where ease of manipulation and adaptability to surface contours is essential.

Long-Lasting

- Shown to provide a stable, long-lasting barrier during healing and integration of Puros[®] Allografts and staged or immediately placed implants.^{1,2}

Retains Natural Tissue Structure

- Supports an aesthetic soft tissue response^{1,2} through facilitation of cell attachment and proliferation and remodeling into vascularized connective tissue.^{1,3}
- Retains the structure and composition of natural tissue due to the proprietary Tutoplast process, leading to optimal performance and handling.^{1,2}



The Unique Tutoplast Process

More than 11 million implants have been sterilized through the Tutoplast Process with zero confirmed incidence of implant associated infection.⁴


The proprietary Tutoplast process assures the highest standard of tissue safety and quality.⁵ The process preserves the valuable collagen matrix and tissue integrity while inactivating pathogens and gently removing unwanted materials, such as cells, antigens and viruses.⁵ The result is safe, biocompatible tissue.




1. Alkaline Treatment




2. Osmotic Treatment



3. Oxidative Treatment



4. Solvent Treatment



5. Irradiation

Take a closer look



Immediate extraction sites, maxillary left



Immediate Tapered Screw-Vent® Implant placement and grafting with Puros Cancellous Particulate Allograft



CopiOs Pericardium Membrane placement



CopiOs Pericardium Membrane in combination with autologous fibrin



Six-month re-entry

1. Kistler S, Bayer G, Kistler F, Am Lech L. Experience with the biological Tutodent membrane in implant practice. *Implantologie Zeitung Journal*. 2004;8(7):47-48.
2. Simsek B, Simsek S. Evaluation of success rates of immediate and delayed implants after tooth extraction. *Chinese Medical Journal*. 2003;116(8):1216-1219.
3. Steigmann M. Pericardium membrane and xenograft particulate grafting materials for horizontal alveolar ridge defects. *Implant Dent*. 2006;15:186-191.
4. Data on file with RTI Surgical, Inc.
5. Tadic D, Epple M. A thorough physicochemical characterization of 14 calcium phosphate-based bone substitution materials in comparison to natural bone. *Biomaterials*. 2004;25:987-994.

Photos ©2012 Dr. Maurice A. Salama, DMD. All rights reserved. Individual results may vary.

For more information, visit ZimVie.com

ZimVie
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Phone: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV1967 REV A 12/23 ©2023 ZimVie. All rights reserved.

Ordering Information

Catalog Number	Description
77776	CopiOs Pericardium Membrane, 15 x 20 mm
77777	CopiOs Pericardium Membrane, 20 x 30 mm
77778	CopiOs Pericardium Membrane, 30 x 40 mm

ZimVie
Restoring Daily Life.™

